C3

INSTRUCTION FOR USE

ENGLISH	2
FRANÇAIS	7
ESPAÑOL	12
DEUTSCH	17
ITALIANO	22

SpineVision® C³® Anterior Cervical Plate System (C³® ACPS) Information

Purpose

The C³® Anterior Cervical Plate System (C³® ACPS) is intended as a temporary internal fixation device used for the correction and stabilization of the cervical spine.

The system is also intended to enhance the development of a solid fusion.

According to the clinical evaluation report, it can be concluded that the benefice-risk ratio is positive.

Description

The C³® ACPS consists of a number of plates and screws of varying lengths and a variety of cover plates with locking set screws. There are also ancillary products. All screw and plate components are machined from implant grade titanium alloy,

Ti-6Al-4V(ELI) per ISO 5832-3 and ASTM F-136.

Indications and Usage

Properly used, the C³® ACPS will aid in the development of a solid spinal fusion. The surgeon must be aware of not only the medical and surgical aspects of the implant, but also of the material limitations of metallic implants. The patient must be instructed in the limitations of the implant, especially with regard to the weight bearing and other body stresses on the device before a solid fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device. The patient must be instructed that device failure could lead to the need for additional surgery to remove failed components.

The C³® ACPS is indicated for trauma, deformity (lordosis, kyphosis and scoliosis), pseudoarthrosis, previously failed cervical spine fusion, tumor, degenerative disc disease (DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies) spondylolisthesis, and spinal stenosis.

The C³® ACPS is indicated for stabilizing the spine from C2 to C7. Unicortical or bicortical fixation from the anterior face of the vertebral bodies may be employed using either the variable angle or fixed angle screws. This product may be employed as an internal fixation device during the time interval required for arthrodesis.

Packaging:

The packaging of each component must be intact on reception. The operating room personnel must thoroughly verify, before use, that all devices are complete and that no component presents any signs of damage. Damaged packagings and/or products must not be used, but must be returned to SpineVision®.

All packaging and labeling materials must be removed prior to sterilization. Only sterile implants and instruments must be used in surgery.

Cleaning and Sterilization:

The C³® ACPS instruments and implants are supplied not sterile. They must be sterilized before use. Surface decontamination and cleaning of reusable instruments must be completed prior to sterilization.

For the initial and any subsequent cleaning, the recommended general protocol for the pre-disinfection, decontamination and cleaning of **the implants and instruments** is:

Disassemble necessary implants and instruments. Articulated instruments must be opened.

In a pre-disinfecting bath, immerse and soak (minimum of 15 minutes) implants and instruments in a quaternary ammonia detergent. Transfer implants and instruments to an ultrasound bath, immerse and soak (15 minutes) implants and instruments under ultrasound action, in a quaternary ammonia detergent. Detergents must be used according to manufacturer recommendations. It is recommended to use deionized or distilled warm (room temperature) water for soaking, cleaning and rinsing.

After soaking, scrub all surfaces with a soft appropriate brush, paying close attention to threads, cannulas, hinges and hard to reach areas.

Rinse immediately and thoroughly after washing.

Inspect each device and repeat the washing in the presence of any residual debris.

Immediately dry products.

Note: Certain cleaning solutions such as those containing bleach or formalin may damage some devices, and should not be used, except in the case of aluminium instruments used on high risk patients (see Note below).

Lubrication specification:

SpineVision® recommends the lubrication of instruments that are mechanical or contain articulating surfaces. The lubricant used for the lubrication process of surgical instruments must be available in spray and is recommended to bear the CE marking according to the 93/42/EEC Directive modified by 2007/47/EEC Directive. The instruments must be lubricated after each pre-disinfection/decontamination and before each steam sterilization.

Prior to use, the C³® ACPS System components should be steam sterilised according to the following parameters:

METHOD	CYCLE	TEMPERATURE	MINIMUM EXPOSURE TIME
Steam	Vacuum	134°C (273°F)	18 minutes

The instruments and implants must be completely dried before use.

Verification:

Devices must always be verified before use: any devices presenting signs of weakness or surface scratches must not be used.

Note:

New requirements related to high risk patients:

In the case of high risk patients, i.e. being suspected of contact with Non-Conventional Transmissible Agents (e.g. prions), decontamination should be performed according to requirements of the World Health Organization. In this case, SpineVision® instruments can be decontaminated with sodium hydroxide, with the strict exception of instruments made out of aluminum, where bleach should be used.

Postoperative mobilization

Careful patient handling for the first two or four months following the operation is very important while the fusion mass matures and becomes able to share load with the implant. Instructing the patient to reduce stress on the implants is an important part of the effort to avoid clinical problems that may accompany failed fixation.

External immobilization is recommended until the maturation of the fusion mass is confirmed radiographically.

Once the fusion mass has healed, the device should be removed to prevent stress-shielding of the fused bone.

Contraindications

- widely disseminated metastatic tumors of adjacent vertebral bodies.
- severe osteoporosis.
- overt infection of the involved vertebral bodies.
- any entity or condition which totally precludes the possibility of successful fusion such as malignant tumour, kidney dialysis, congenital anomalies, or osteopenia.
- other relative contraindications include, but are not limited to, obesity, certain degenerative diseases, and foreign body sensitivity (known or suspected allergy or intolerance to metals). The patient's activity level, mental condition, or occupation may be factors that affect the surgical outcome.
- Other conditions, including alcoholism and drug abuse may also place excessive stresses on the device.
- rapidly evolving joint diseases, bone absorption
- do not use this device in the presence of any neural or vascular deficits or compromising pathology that may be further injured by device application.
- all patients not agreeing to comply with postoperative instructions.

Warnings

- correct selection of the implant is extremely important. The chance of adequate fixation is increased by the selection of the proper size of the implant. Proper selection of the implant can minimize, but not eliminate risks. No internal fixation device can withstand activity levels equal to those withstood by normal healthy bone. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.
- mixing metals can cause corrosion. Titanium and stainless steel in contact will accelerate the corrosion process in the stainless steel. The presence of corrosion may accelerate the fatigue failure of implanted devices. The amount of metallic compounds released into the body will be increased. Therefore, internal fixation devices which may come into contact with other metal devices must be made from similar or compatible materials.
- the C³® ACPS components have been designed for use with this system. Do not use the C³® ACPS with any components not specifically recommended.
- the C³® ACPS is not approved for screw attachment or fixation to the posterior elements (pedicles) of cervical, thoracic or lumbar spine.
- implants may fail when subjected to the increased loading associated with delayed union or nonunion. Internal fixation devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur

or is delayed, the device may break due to fatigue. Notches, scratches, or over-bending of the implant during the course of surgery may contribute to early failure. Patients should be informed of the risks of implant failure.

- patient selection : the following factors may be of importance when selecting patients for internal fixation devices :
 - foreign body sensitivity
 - certain degenerative diseases
 - senility, mental illness, alcoholism, or drug abuse
 - obesity
 - patient's occupation or activity level

General Precautions

- **SURGICAL IMPLANTS MUST NEVER BE REUSED;**
An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage.
- **correct handling of the implants is extremely important.** Care must be taken when contouring the plate components of the C³®ACPS. Alterations may produce defects in surface finish and may create internal stresses that could become the focal point for eventual breakage of the implant.
- **adequately instruct the patient in the appropriate postoperative care.** The patient's ability and willingness to follow instructions are one of the most important aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal fixation devices by loosening, bending, or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may occur if excessive loading is placed on the implant. An active, debilitated, or demented patient who is unable to use orthotic weight bearing devices may be particularly at risk.
- **removal of the implant after healing.** Particularly in young active patients, implants may loosen, fracture, corrode, migrate, increase the risk of infection, cause pain, or stress shield bone – even after normal healing. The surgeon should consider the risks and benefits when deciding whether or not to remove an implant. Careful postoperative management to avoid re-fracture should follow implant removal. If the patient is older and has a low activity level, the surgeon may elect not to remove the implant in order to eliminate the risks of another surgery.
- **proper placement of the implant.** The implant should be placed as close as possible to the midline of the spine.
- **proper fastening of screws.** It is important that the surgeon follows the recommendations for screw insertion and tightening. Insufficient fastening of the implant components may result in premature device failure. Over-tightening screws may cause the thread to strip in the bone, and may result in loosening and potential implant failure.

Possible Adverse Affects

- delayed union
- non union
- bending, fracture, loosening, or migration of implant

- metal sensitivity or allergic foreign-body reaction
- infection
- stress shielding resulting in decreased bone density
- pain, discomfort, or other abnormal sensations due to the presence of the device
- erosion of blood vessels due to the proximity of the device. Such erosion may lead to hemorrhage which could result in death.
- bursitis
- dural leak
- paralysis
- death

Product complaints

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identify, durability, reliability, safety, effectiveness and/or performance, should notify SpineVision®. Further, if any of the implanted C³® ACPS components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately.

If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filling a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

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Email: corp.quality@spinevision.com

IMPORTANT INFORMATION CONCERNING C³® ACPS System INSTRUMENTATION

Caution: Federal Law (USA) Restricts This Device To Sale By Or On The Order Of A Physician.

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SPINEVISION Inc.