Flex +2

INSTRUCTION FOR USE

ENGLISH	2
FRANÇAIS	7
ITALIANO	12
ESPAÑOL	17
PORTUGUÊS	22
DEUTSCH	28

SpineVision®: Flex⁺²

Introduction

Flex⁺² implants are designed for the correction and surgical stabilization of the spinal column only in combination with pedicle screw systems manufactured and marketed by SpineVision : titanium P.L.U.S.®, titanium X-PLUS®, LUMIS™ or U.L.I.S.™. The rods have been designed to allow the surgeon to adapt the stabilization type to the degree of degeneration.

According to the clinical evaluation report, it can be concluded that the benefice-risk ratio is positive.

Description

The Flex⁺² range comprises different types of rods: dynamic stabilization rods (F1-R1609T), and hybrid stabilization rods allowing the surgeon to associate an arthrodesis with a dynamic stabilization thanks to the flexible part of the rod on one level (F1-RH1609T, F1-RHL1609T, F1-RH1609CT) or two contiguous levels (F1-RH2609T).

The metallic components are made of titanium alloy (Ti-6Al-4V) compliant with standard ASTM F136 (ISO 5832-3), and the parts used for dynamic stabilization are overmolded with biocompatible Urethane PolyCarbonate (USP class VI-compliant).

WARNING: the Flex⁺² rods are not to be used with screws or instruments from another manufacturer, or with pedicle systems not referenced in the instructions.

Indication, use and contraindications

Indications

When properly used on a patient whose skeleton is mature, **any flexible part of the Flex⁺²** rod is intended to stabilize vertebral segments over one or two adjacent stages from T10 to L5 in the following spinal pathologies :

- Moderate lumbar discopathy at Pfirrmann stages III or IV
- Early degenerative spondylolisthesis (grade I) and hyper-mobile segment without extensive facet joint osteoarthritis
- Lumbar stenosis by bulging of discs and of flaval ligaments, associated with a recalibration
- Massive herniated disc in a young patient with anticipated discopathy

Any rigid part of the Flex⁺² rod is intended to fuse several vertebral segments from T11 to S1, directly adjacent to the dynamically stabilized level.-in the following spinal pathologies :

- Degenerative discopathies
- Facet degeneration
- Degenerative spondylolisthesis
- Lumbar stenosis
- Pseudarthrosis following a previous instrumented surgery

Contraindications

Contraindications include but are not limited to :

- All cases not described in the indications.
- Treatment of thoracic spine above T10 and of cervical spine.
- Unilateral application of the Flex⁺².
- Use of two rods with different references on the same level

- Use of P.L.U.S.[®] titanium screws, X-PLUS[®] titanium screws, LUMIS[™] or U.L.I.S.[™] screws with a diameter below 5.5 mm at the flexible level,
- Any previous surgery at the flexible level(s) or previous instrumented lumbar surgery at the rigid level(s) except if pseudarthrosis can be proved,
- Pathological obesity (Body mass index, BMI > 35).
- Pregnancy.
- Concomitant performance of facetectomy and/or laminectomy at the level(s) instrumented with the flexible part(s)
- Known or suspected allergy, metal and/or urethane polycarbonate intolerance.
- Localized infection of the surgical site, signs of inflammation, fever, systemic or local active infection.
- All patients with insufficient tissue coverage of the surgical site.
- Rapidly developing arthropathy, acceleration of bone resorption, osteopenia / osteoporosis or osteomalacia highlighted by X-rays or fracture antecedents. The presence of osteoporosis is a relative contraindication, as it can limit the expected correction gain and the stability of the mechanical fastening.
- Smokers whose bone density is below a T-Score of -1.5.
- Arthrosis at the flexible level(s)
- Degenerative scoliosis > 10° at affected level(s).
- Spondylolisthesis higher than grade I and spondylolisthesis by spondylolysis at the flexible level(s)
- Any other medical or surgical affection liable to compromise the success of the instrumented operation, such as, extensive deformation, tumor, vertebral fractures or serious congenital anomalies, acceleration of the erythrocyte sedimentation speed not accounted for by another disease, high leukocyte count or tendency to leukopenia.
- All cases requiring the association of different metals.
- Mental illness.
- Any patient who does not accept to follow the postoperative instructions.

The contraindications for these devices are similar to those for other instrumentations by vertebral spinal rods. This spinal instrumentation device is not designed or sold for applications other than those mentioned.

Warning: Use of the bender on the flexible part(s) of the Flex⁺² implant is prohibited to avoid damaging the overmolded part. For hybrid rods, the flexible part must always be implanted at the cranial end of the assembly.

The surgeon must have a sound knowledge, not just of the medical and surgical aspects of the implant, but also of the material limits of implants in urethane polycarbonate.

The patient must be informed of the limits of the implant, particularly in terms of loads and other bodily stresses exerted on the device prior to forming of a solid arthrodesis with fused stages.

The patient must also be aware that failure to comply with the postoperative instructions could lead to failure of the device, possibly requiring a new operation to remove it.

This device can only be implanted by a surgeon who is thoroughly familiar with it and with its applications, and with the instruments and the surgical technique required.

Patient selection

- Particular attention must be paid to the patient selection criteria, to the proper setting of the implant, and to postoperative care, in order to minimize the stresses to which the implant is subjected. In particular, it is important to take into consideration the patient's sagittal equilibrium so as to take his pelvic equilibrium parameters into account and avoid excessive shear stresses on the flexible part(s).
- Only the patients meeting the criteria described in the indications can be selected.
- Patients who meet the criteria described in the contraindications must not be selected.

Possible adverse effects

The following adverse events are inherent to any surgical operation, or associated with instrumented intervention on the spine, including in particular:

- Coming apart, warping, slipping and/or rupture of any or all of the components or instruments during surgery.

- Reaction to the implant as a foreign body, e.g. mass, autoimmune disease, metallosis and/or healing complications.
- Decrease in curvature or correction of the spinal column or in the height of the disc.
- Infection.
- Fracture of the vertebral body due to load transfer above and below the instrumented area.
- Absence of bone setting (pseudarthrosis) at the rigid part.
- Pressure on the skin or other soft tissues from component parts possibly causing skin penetration, irritation, fibrosis, necrosis, and/or pain.
- Device breakage, possibly triggering pain
- Loosening and migration of any or all of the components
- Alteration of neurological functions, occurrence of radiculopathies, dural lesions and/or pain,neurovascular deficit, including paralysis or other serious lesions. Effusion of cerebrospinal fluid.
- Digestive, urological and/or genital disorders, including sterility, impotence.
- Hemorrhage and/or hematoma.
- Disk inflammation, arachnoiditis and/or other types of inflammation.
- Deep venous thrombosis, thrombophlebitis and/or pulmonary embolism.
- Complications at bone graft donor site.
- Inability to resume usual day-to-day activities.
- Lesion of soft tissues.
- Death

NOTE: Certain potential adverse events may require an additional surgical operation.

Storage and packaging

The packaging of each of the components must be intact on delivery. A “tamper-proof” safety label testifies to the integrity of the packaging. If the word OPENED is visible, then the packaging must be considered as damaged, and is not to be used.

The sterility use-by date must be complied with. Make sure that the implants are not damaged when opening the packaging. Damaged packages or products must not be used. If in doubt, please send them back to SpineVision.

Cleaning and sterilization

Flex⁺² implants are delivered sterile (Ethylene Oxide [EO marking on the packaging]). They must not be resterilized before use.

Flex⁺² implants must only be used with the titanium P.L.U.S.[®], titanium X-PLUS[®], LUMIS[™] or U.L.I.S.[™] ancillary equipment. The titanium P.L.U.S.[®], titanium X-PLUS[®], LUMIS[™] or U.L.I.S.[™] implants and instruments are provided clean but not sterile. All packaging and labeling equipment must be removed before sterilization. Decontamination and cleaning must be performed before sterilization.

For initial or future cleaning, the general protocol recommended for pre-disinfection, decontamination and cleaning of the **instruments** is as follows:

Dismantle the necessary instruments. Hinged instruments must be opened.

In a pre-disinfection tank, immerse the instruments and allow them to soak (15 minutes) in a detergent solution of quaternary ammonium compound. Transfer the instruments into an ultrasonic tank, immerse the equipment and allow soaking (15 minutes). The detergent must be used in accordance with the manufacturer’s recommendations. It is recommended to use demineralized or distilled water at ambient temperature for immersion, cleaning and rinsing.

After immersion, brush all surfaces using an appropriate soft brush, taking particular care with the threads, cannulae, hinges and hard-to-reach areas.

Immediately and thoroughly rinse after cleaning.

Check each instrument and clean again if there are any residues remaining.

Dry the instruments immediately after rinsing.

Note: Certain cleaning solutions, such as those containing bleach or formaldehyde, may damage certain components and must not be used, except for aluminum instruments used on high-risk patients (see note below).

Lubrication specifications

SpineVision recommend lubricating instruments containing gears or hinged surfaces. The lubricant used for surgical instruments must be available in spray form, and it is recommended that it should bear the EC marking as per directive 93/42/CEE. The instruments must be lubricated after each pre-disinfection / decontamination and before each steam sterilization.

Sterilization

Before use, titanium P.L.U.S.[®], titanium X-PLUS[®], LUMIS[™] or U.L.I.S.[™] instruments and implants must be sterilized with steam in accordance with the following parameters:

Method	cycle	temperature	Minimum exposure time
Steam	Vacuum	134° (273°F)	18 minutes

The instruments must be thoroughly dried before use.

Verification

The devices must always be checked before use: any item showing signs of fatigue or surface cracks must not be used.

Note: In the case of a high-risk patient, i.e. suspected contact with Unconventional Transmissible Agents (e.g. prions), the decontamination must be carried out in accordance with the requirements defined by the World Health Organization (WHO). In this particular case, the Spinevision instrumentation can be decontaminated with sodium hydroxide, with the strict exception of instruments made from aluminum, for which bleach must be used.

Precautions

Pedicle screw spinal systems must be implanted by an experienced spinal surgeon specifically trained in their use, as the technique required for their implantation is demanding, and exposes the patient to a serious risk of lesion.

Preoperative precautions:

- Implants are to be handled and stored with the utmost care. They must not be scratched or damaged. The implants and instruments must be protected while in storage, particularly against corrosive environments.
- The surgeon must be familiar with all the components before using the equipment, and must personally make sure that all the implants and instruments required are available before the start of the operation.
- The type of construction to be performed in each case must be determined before the start of the operation. An appropriate range of implants of various sizes, including smaller and larger than those initially provided for, must be available at the time of the operation,
- All components must be sterile before they can be used. Extra sterile components must be available, should they be needed unexpectedly.

Peroperative precautions:

- The surgeon must strictly follow the instructions for use of spinal instrumentation.
- The surgeon must be extremely careful to avoid any lesion to the spinal cord or to the nerve roots. This warning is of particular importance when inserting the screws. Any nervous lesion may induce a neurological dysfunction.
- A rupture, slippage or improper use of the instruments or implant components may cause injury to the patient or to members of the surgical team.
- The rigid parts of the Flex⁺² rods must not be repeatedly or excessively bent beyond what is strictly necessary.
- The slightest stress on the Urethane PolyCarbonate part of Flex⁺² rods is prohibited, whether during bending of the rod or while it is being implanted.
- In particular, make sure no pre-stressing is applied to the Urethane PolyCarbonate part of the Flex⁺² rods between the screws.
- The screw tightening components must not rest on the Urethane PolyCarbonate part, as this might damage the part.
- Check very carefully that the surfaces of the implant are neither scratched nor damaged before implantation, and in particular, check that the jointing between the titanium rings and the Urethane PolyCarbonate part is not damaged.

- If the rods are to be cut to a certain length, the section must be a flat non-abrasive surface perpendicular to the axis.
Cut the rods outside the surgical field.
- Do not use screws with unsuitable dimensions (length, diameter), as this might damage the nerve roots, or cause hemorrhages and/or tearing.
- Before closing, all locking screws must be tightened as per the instructions.

Postoperative precautions:

- The postoperative recommendations given by the surgeon to the patient, and compliance with these recommendations by the patient, are extremely important.
- The patient must be informed of the limitations in the use of the device. He must be aware that complications such as deformation, coming apart or rupture of the device may occur if the spinal column is subjected to premature and/or excessive loads.
- The patient must be aware that motion of the spinal column will be limited at the instrumented level, and will have to be trained in making up for this permanent physical limitation in bodily motions.
- A persistent absence of bone setting at the fusion level which cannot be immobilized exerts repeated and excessive stresses on the implant. Due to the material fatigue phenomenon, these stresses may cause the device to come apart, be distorted or ruptured. It is important to mobilize the fusion area and to confirm this setting through an X-ray examination. In case of abnormal absence of setting or of coming apart, distortion and/or rupture, the device is to be removed immediately, before it causes a serious lesion.

All explanted devices must be treated so as to make it impossible for them to be re-used during another surgical operation.

SURGICAL IMPLANTS MUST NEVER BE RE-USED. Once removed from a patient's body, an implant must never be re-implanted. Even if a device seems intact, it may present small defects and internal fatigue traces liable to cause premature rupture.

Complaints about a product

Any health professional (e.g. prescribing physician or user of the systems described in this document) who wishes to lodge a complaint, or who is not satisfied with the qualities, identification, durability, reliability, safety of use, efficiency and/or performances of the product, is invited to inform SpineVision . Moreover, SpineVision must immediately be informed if one of the implanted components has not performed as specified, or does not work as expected.

SpineVision must be informed as soon as possible by phone, fax, mail or e-mail if any one of its products may have been responsible for the death of, or serious injuries to, a patient, or may have been a contributing factor. If you lodge a complaint, please give the name and the number(s) of the component(s) involved, the batch number(s), your name and address, the nature of your complaint, and state whether you are requesting a written report from SpineVision.

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