## SUBSTITUT OSSEUX SYNTHÉTIQUE FR

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## SYNTHETIC BONE SUBSTITUTE EN

Orthopaedic application - Instructions for use

DESCRIPTION

The Synthetic Bone Substitute is a bone void filler. This product is a microporous and macroporous biphasic calcium phosphate ceramic consisting of 60% Hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP).

The Synthetic Bone Substitute is presented in a porous form required for the biological exchanges particularly for bone ingrowth.

The Synthetic Bone Substitute may be used with physiological saline, patient's own serum, whole blood, or bone marrow aspirate (BMA).

The Synthetic Bone Substitute is available in various shapes packaged in vials or syringes to be easily adapted to several bone defect designs. This product is provided sterile for single patient use.

The Synthetic Bone Substitute is progressively replaced by new bone according to the remodeling process.

INTENDED USE

The Synthetic Bone Substitute is intended for use to reconstruct bony voids or bone gaps of the skeletal system (e.g. extremities, spine and pelvis).

The performances of the Synthetic Bone Substitute are the filling of bone defects and the bony ingrowth from local osseous tissue onto the surface of the product (osteoconduction process). The total resorption of the product depends on numerous factors such as its size and volume, the location of the defect, the surgical technique, the health status of the patient. The bone reconstruction time is always inferior to the total resorption time.

SURGICAL INDICATIONS

The Synthetic Bone Substitute is indicated for bony voids surgically created or bone gaps osseous defects created from traumatic injury to the bone.

The Synthetic Bone Substitute can be used with autograft as a bone void filler extender.

The target population is adults, excluding pregnant women. The Synthetic Bone Substitute has not been tested on pregnant women and neither for use in pediatric population.

CONTRAINDICATIONS

The Synthetic Bone Substitute is contraindicated where the device is intended as structural support in the skeletal system.

Conditions representing contraindications include:

- Ostomyelitis
- Infection in necrotic surgical sites
- Degenerative bone disease
- Intra-articular implantations
- Opening meninx
- Long-term steroidal therapy or treatment acting on the calcium or phosphorus metabolism.

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Some of adverse effects may be treated by drugs or require reintervention.

STERILISATION AND CONSERVATION

The Synthetic Bone Substitute is provided sterile and is sterile unless the inner packaging has been opened or damaged. The Synthetic Bone Substitute must not be resterilized. The sterilization process of this product using radiations was validated during its conception. We do not guarantee the sterility of the product if a new sterilization is performed.

The Synthetic Bone Substitute is for single patient use only and should never be reused. The non-respect of the single use can conduct to septic issues.

The Synthetic Bone Substitute must be stored in its original packaging and kept at a temperature between 15°C and 30°C (59°F and 86°F).

Any product or container which has been in contact with the patient or is contaminated, should be treated as clinical waste and disposed of in accordance with the local regulations for the disposal of clinical waste.

**WARNING AND PRECAUTIONS**
**DO NOT USE IF THE PACKAGE HAS BEEN OPENED OR DAMAGED OR IF THE EXPIRATION DATE HAS BEEN EXCEEDED.**

Do not overfill defects.

The maximum amount per site of implantation is 50 cm³.

The device should be secured to prevent motion and migration, used in areas where the graft can be adequately contained.

Do not leave defect open.

As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure.

This device does not render after implantation until healing. This device is a void filler used in cavity with respected corticals or associated with an osteosynthesis device.

Do not compromise blood supply to the defect area.

Maximum contact between the product and the recipient bone must be established.

This device has to be in direct contact of at least 3 surfaces.

Adequate post-operative immobilization is necessary to prevent movement that may lead to soft tissue ingrowth.

The Synthetic Bone Substitute should be manipulated with precaution in order to preserve the porous structure, it must not be cut or shaped before implantation.

The Synthetic Bone Substitute is intended for use by physician familiar with bone void filling and rigid fixation techniques.

The Synthetic Bone Substitute is to be handled or implanted by trained qualified physicians having read this instruction for use.

**Caution:** Federal law in USA and Canada restricts this device to sale by or on the order of a physician.

PRÉPARATION

- Use care to avoid destruction of porous structure.
- Prior to this, product must be hydrated without excess with physiological saline.
- This product may be soaked in patient's serum, whole blood, or bone marrow aspirate (BMA) until graft is completely impregnated.
- The proportion of fluid should be 1 volume fluid to 2 volumes of product

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