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### DESCRIPTION

The Injectable Bone Substitute is a bone void filler. This product is a micro porous and macroporous bi-phase calcium phosphate ceramic consisting of 60% Hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP) particles associated with an excipient vehicle of pharmaceutical grade quality (hydroxypropylmethylcellulose in an aqueous solution). The soluble carrier acts as a spacer and a binder of the particles. This device provides an environment for bone ingrowth. The Injectable Bone Substitute is provided sterile in a syringe.

The Injectable Bone Substitute is progressively replaced by new bone according to the remodeling process.

**INTENDED USE**
The Injectable Bone Substitute is intended for use to reconstruct bony voids or bone gaps of the skeletal system (e.g. extremities, spine and pelvis).

The performances of the Injectable Bone Substitute are the filling of bone defects and the bony ingrowth from local osseous tissue onto the surface of the product (osteoconduction process).

The total resorption of the product depends on numerous factors such as its size and volume, the location of the defect, the surgical technique, the health status of the patient. The bone reconstruction time is always inferior to the total resorption time.

### SURGICAL INDICATIONS

The Injectable Bone Substitute is indicated for bony voids surgically created or bone gaps osseous defects created from traumatic injury to the bone.

The Injectable Bone Substitute can be used with autograft as a bone void filler extendor.

The target population is adults, excluding pregnant women. The Injectable Bone Substitute has not been tested on pregnant women and neither for use in pediatric population.

### CONTRAINDICATIONS

The Injectable Bone Substitute is contraindicated where the device is intended as structural support in the skeletal system.

Conditions representing contraindications include:

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### Local :

- Trophic trouble
- Un-treated infection
- Opening meningis
- Irradiated patient
- Granuloma or non-cured cyst
- Cellulitis
- Osteomyelitis aiguë ou chronique
- Intra-articular implantation
- Implantation in necrotic surgical sites

### ADVERSE-EFFECTS

- Post-operative discomfort
- Infection
- Inflamation
- Nécrose
- Infection
- Réaction à corps étranger
- Douleur
- Flèvre
- Réaction inflammatoire
- Adverse effects on blood coagulation and blood components
- Thrombosis
- Necrosis
- Risk of allergy
- Non-union/Pseudarthrose due to displacement product or fibrous encapsulation
- Pain
- Fever
- Degradation of organ

Some of adverse effects may be treated by drugs or require reintervention.

### STERILIZATION AND CONSERVATION

The Injectable Bone Substitute is supplied sterile in a syringe. This device is sterilized by steam. This device should not be resterilized. In all operations involving this device, appropriate aseptic technique must be observed.

The Injectable Bone Substitute is for single use only and should never be resterilized. The sterilization process of this product was validated during its conception.

We do not guarantee the sterility of the product if a new sterilization is performed on it. The non-respect of the single use can conduct to septic issues.

The Injectable Bone Substitute is for single patient use only and should never be reused. The non-respect of the single patient use can conduct to septic issues.

The Injectable Bone Substitute must be stored in its original packaging and kept at a temperature between 15°C and 30°C (59°F and 86°F).

Any product or container which has been in contact with the patient or is contaminated, should be treated as clinical waste and disposed of in accordance with the local regulations for the disposal of clinical waste.

**WARNING AND PRECAUTIONS**
DO NOT USE IF PACKAGING IS OPENED OR DAMAGED OR IF EXPIRATION DATE HAS BEEN EXCEEDED.

The maximum amount per site of implantation is 40 cm<sup>3</sup>.

The injectable Bone Substitute should be secured to prevent motion and migration, used in areas where the graft can be adequately contained.

Do not leave defect open.

As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure.

The Injectable Bone Substitute does not harden after implantation until healing. This device is a void filler used in closed cavity with respected cortical or associated with an osteosynthesis device. Maximum contact between the product and the recipient bone must be established. This device has to be in direct contact of at least 3 surfaces. Do not compromise blood supply to the defect area.

Adequate post-operative immobilization is necessary to prevent movement that may lead to soft tissue ingrowth.

Until healing is complete this device may not withstand weight bearing. Post-surgery discharge is necessary.

The implantation in a revision surgical site containing non-resorbable fragments of material (e.g. polyethylene ligament waste, carbon fibers) is not recommended.

The delivery ancillary provided for reference less than 3ml, is made of glass. A special care must therefore be taken: the syringe cannot be used for other ancillary purposes.

To avoid the extraction of the material, the filling of fenestrated cage is not recommended. Any excessive washing can cause degradation of the material.

The Injectable Bone Substitute is to be handled or implanted by trained qualified physicians having read this instruction for use.

**Caution:** Federal law in USA and Canada restricts this device to sale by or on the order of a physician

**PREPARATION**
Because of its particular syringe, the product requires a preparation advice explained in leaflet PRCR19.03.

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