# Spacevision

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## **SpineVision® Spacevision™ ACIF: Anterior Cervical Interbody Fusion Cage**

### **Purpose**

The Spacevision™ ACIF Anterior Cervical Interbody Fusion Cage is intended as an internal spacer between the vertebral bodies of the cervical spine used in surgical procedures to correct and stabilize the cervical spine.

The system is also intended to enhance the development of a solid fusion.

According to the clinical evaluation report, it can be concluded that the benefice-risk ratio is positive.

### **Description**

The Spacevision™ Set consists of cages of varying heights, widths and depths and ancillary products for placement of the cages. All cages are machined from implant grade PEEK-OPTIMA® LT1.

### **Indications and Usage**

Properly used, the Spacevision™ ACIF is intended to aid in the development of a spinal fusion. The surgeon must be aware of not only the medical and surgical aspects of the implant, but also of the material limitations of PEEK-OPTIMA® (Polyetheretherketone) implants. The patient must be instructed in the limitations of the implant (see Warnings and Precautions), especially with regard to the weight bearing and other body stresses on the device before a solid fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device. The patient must be instructed that device failure could lead to the need for additional surgery to remove failed components.

The Spacevision™ ACIF cage is indicated for degenerative disc disease (DDD is defined as neck pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies), spinal stenosis and previously failed cervical spine fusion (pseudoarthrosis).

### **Packaging**

Packages for each of the components should be intact upon receipt. Care should be taken not to damage the implants when opening the packaging. Damaged packages or products must not be used, but should be returned to SpineVision®.

### **Cleaning and Sterilization**

The Spacevision™ ACIF instruments and implants are supplied not sterile. They must be sterilized before use. All packaging and labeling materials must be removed prior to sterilization. Decontamination and cleaning must be completed prior to sterilization.

**For the initial and any subsequent cleaning**, the recommended general protocol for the pre-disinfection, decontamination and cleaning of the **implants and instruments** is:

Disassemble necessary instruments. Articulated instruments must be opened.

In a pre-disinfecting bath, immerse and soak (minimum of 15 minutes) implants and instruments in a quaternary ammonia detergent. Transfer instruments in an ultrasounds bath, immerse and soak (15 minutes) implants and instruments under ultrasound action, in a quaternary ammonia detergent. Detergents must be used according to

manufacturer recommendations. It is recommended to use deionized or distilled warm (room temperature) water for soaking, cleaning and rinsing.

After soaking, scrub all surfaces with a soft appropriate brush, paying close attention to threads, cannulas, hinges and hard to reach areas.

Rinse immediately and thoroughly after washing.

Inspect each device and repeat the washing in the presence of any residual debris.

Immediately dry products.

**Note:** Certain cleaning solutions such as those containing bleach or formalin may damage some devices, and should not be used, except in the case of aluminium instruments used on high risk patients (see Note below).

**Lubrication specification:**

SpineVision® recommends the lubrication of instruments that are mechanical or contain articulating surfaces. The lubricant used for the lubrication process of surgical instruments must be available in spray and is recommended to bear the CE marking according to the Directive 93/42/EEC. The instruments must be lubricated after each pre-disinfection/decontamination and before each steam sterilization.

**Sterilization:** Prior to use, the Spacevision™ ACIF instruments and implants should be steam sterilized according to the following parameters:

METHOD	CYCLE	TEMPERATURE	MINIMUM OF EXPOSURE TIME
steam	vacuum	134°C (273°F)	18 minutes

The instruments and implants must be completely dried before use.

Check for evidence of deterioration on the implants after sterilization (return to the manufacturer or distributor any implants showing deterioration).

**IMPORTANT NOTE:** As the PEEK Optima® LT1 glass transition temperature is 145°C (293°F), it is important not to use a higher sterilization temperature than 134°C (273°F). Otherwise the implant could suffer deformation.

To minimize the risk of deformation due to too high temperatures, do not place any load on top of the implants during cleaning and sterilization.

**Note:**

New requirements related to high risk patients:

**In the case of high risk patients**, i.e. being suspected of contact with Non Conventional Transmissible Agents (e.g. prions), decontamination should be performed according to requirements of the World Health Organization. In this case, SpineVision® instruments can be decontaminated with soda, with the strict exception of instruments made out of aluminium where bleach should be used.

Careful patient monitoring for the first two to four months following the operation is very important while the fusion mass matures and becomes able to share load with the implant. The patient must be instructed to reduce stress on the implants in order to avoid clinical problems that may accompany failed fixation.

External immobilization is recommended until the maturation of the fusion mass is confirmed radiographically.

**Contraindications**

- severe osteoporosis, calcium metabolism disorder
- pregnancy
- infection
- other relative contraindications include, but are not limited to, obesity, certain degenerative diseases and foreign body sensitivity. The patient's activity level, mental condition, **or occupation** may be factors that affect the surgical outcome.

- other conditions, including alcoholism and drug abuse may also place excessive stresses on the device.
- do not use this device in the presence of any neural or vascular deficit or compromising pathology that may be further injured by device application.

### **Warnings**

- the use of Spacevision™ ACIF cage is to be done by a trained orthopedic surgeon or a trained neurosurgeon familiar with these recommendations as well as with the operating techniques, relating to these implants and using dedicated adequate instrumentation.
- correct selection of the implant is extremely important. The chance of adequate stabilization is increased by the selection of the proper size of the implant regarding the anatomy and conditions of the patient. Proper selection of the implant can minimize, but not eliminate risks.
- no internal stabilization device can withstand activity levels equal to those withstood by normal healthy bone. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.
- the Spacevision™ ACIF cage components have been designed for use with this system. Do not use the Spacevision™ ACIF cage with any components not specifically recommended, or with any components from another manufacturer.
- the Spacevision™ ACIF cage must not be inserted from the posterior side of the cervical spine nor is it for use or lumbar spine. in the thoracic
- implants may fail when subjected to the increased loading associated with delayed union or non-union. Internal stabilization devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur or is delayed, the device may break due to fatigue. Patients should be informed of the risks of implant failure.
- patient selection: the following factors may be of importance when selecting patients for internal stabilization devices :
  - foreign body sensitivity
  - certain degenerative diseases
  - senility, mental illness, alcoholism, or drug abuse
  - obesity
  - patient's occupation or activity level
  - smoking

### **General Precautions**

- **surgical implants must never be reused.**
- **correct handling of the implants is extremely important.** The implants must be transported in a way not to incur any damage. The implants must not be scratched and must not collide during handling. Care must be taken when picking up the Spacevision™ ACIF cage, filling it with bone graft or bone substitute and inserting it between the

vertebral bodies. Alterations may produce defects in surface finish and may become the focal point for eventual breakage of the implant or for decreased efficacy in withstanding post-operative migration (e.g. when teeth of implant are damaged).

- **adequately instruct the patient in the appropriate postoperative care.** The patient's ability and willingness to follow instructions are one of the most important aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal stabilization devices by migrating, loosening or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may occur if excessive loading is placed on the implant. An active, debilitated, or demented patient who is unable to use orthotic weight bearing devices may be particularly at risk.
  
- **proper placement of the implant.** The implant should be placed on the peripheral cortical bone of the endplates to minimize the risk of post-operative subsidence.
  
- **filling of the implant with bone graft or bone graft substitute.** The Spacevision™ ACIF cage can be filled with bone graft or a bone graft substitute to enhance the bone fusion. The choice of graft will always be the surgeon's responsibility.
  
- **in each case,** the use of an anterior cervical plate, e.g. the SpineVision® C<sup>3</sup>™ ACPS, is recommended in combination with the Spacevision™ ACIF cage.

#### **Possible Adverse Affects**

- delayed union
- non-union
- fracture, loosening, or migration of implant
- infection
- subsidence of implant into adjacent vertebral bodies resulting in loss of height and possible impingement of neural structures
- pain, discomfort, or other abnormal sensations due to the presence of the device
- cage migration that may result in paralysis or blood vessel erosion due to the proximity of the device. Such erosion may lead to hemorrhage which could result in death.
- foreign body reaction due to the presence of implants, such as a mass, auto-immune disease and/or impaired healing
- arrest of growth of fused segments of the spine
- deep vein thrombosis, thrombophlebitis and/or pulmonary embolism
- incapacity to resume the activities of normal everyday life
- death

#### **Product complaints**

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identification, durability, reliability, safety, effectiveness and/or performance, should notify SpineVision®. Further, if any of the implanted Spacevision™ ACIF cages components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately.

If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

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