Spacevision

OLIF PLIF TLIF

INSTRUCTION FOR USE

UNITED STATES OF AMERICA.....	2
Outside U.S.A. - ENGLISH VERSION.....	6
Outside U.S.A. - FRENCH VERSION.....	12
Outside U.S.A. - GERMAN VERSION.....	16
Outside U.S.A. - ITALIAN VERSION.....	21
Outside U.S.A. - SPANISH VERSION.....	25

UNITED STATES OF AMERICA

SpineVision® SpaceVision® PLIF: Posterior Lumbar Interbody Fusion Cage

SpineVision® SpaceVision® TLIF: Transforaminal Lumbar Interbody Fusion Cage

SpineVision® SpaceVision® OLIF: Oblique Lumbar Interbody Fusion Cage

Purpose

The SpaceVision® PLIF (Posterior Lumbar Interbody Fusion), the SpaceVision® TLIF (Transforaminal Lumbar Interbody Fusion) and SpaceVision® OLIF (Oblique Lumbar Interbody Fusion) systems are intended as intervertebral spacers for the lumbar spine (L2-S1). They are intended to be implanted via a posterior approach to give a primary stability of the anterior column (with or without restoring disc height) until fusion occurs.

The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF spacers must be combined with an appropriate posterior lumbar fixation device, e.g. the SpineVision® UNI-Thread™, PLUS®, ULIS™ or LUMIS™.

Description

The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF systems consist of cages of various widths and heights, which can be inserted between two vertebral bodies of L2-S1 segment to give support and correction during lumbar interbody fusion surgeries. The hollow geometry of the implants allows them to be packed with autogenous bone graft.

For the SpaceVision® PLIF system, two cages must be used per segment for fusion in order to stabilize the segment concerned.

For the SpaceVision® TLIF and SpaceVision® OLIF systems, one cage must be used per segment for fusion in order to stabilize the segment concerned.

The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF must be filled with autogenous bone graft or a bone graft substitute to enhance the bone fusion. The choice of graft will always be the surgeon's responsibility. All cages are machined from implant grade PEEK (Polyetheretherketone compliant with ASTM F2026).

Indications

The Spacevision® PLIF (Posterior Lumbar Interbody Fusion), the Spacevision® TLIF (Transforaminal Lumbar Interbody Fusion) and Spacevision® OLIF (Oblique Lumbar Interbody Fusion) systems are intervertebral body fusion devices indicated for use with autogenous bone graft in skeletally mature patients with Degenerative Disc Disease (DDD) at one or two continuous levels from L2-S1. DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies. These DDD patients may also have up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). Patients should have received at least 6 months of non-operative treatment prior to treatment with Spacevision® PLIF, Spacevision® TLIF or Spacevision® OLIF systems. These devices are to be filled with autogenous bone graft material. These devices can be implanted via posterior or transforaminal approaches. Spacevision® PLIF, Spacevision® TLIF or Spacevision® OLIF systems must be used in combination with supplemental internal spinal fixation which has been cleared by the FDA for use in the lumbar spine.

Note for US audience only: The supplemental internal spinal fixation system used has to be cleared by FDA for the corresponding indications.

Contraindications

This device is not intended for thoracic or cervical spine use.

Contraindications include, but are not limited to:

- infection, local to the operative site;
- signs of local inflammation;
- fever or leukocytosis;
- morbid obesity;
- pregnancy;

- pediatric cases, or patient still having general skeletal growth;
- spondylolisthesis unable to be reduced to Grade I;
- suspected or documented allergy or intolerance to composite materials;
- any case where the implant components selected for use would be too large or too small to achieve a successful result.
- any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- prior fusion at the level to be treated.
- any case not needing a bone graft or fusion;
- any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis;
- any other condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors or congenital abnormalities, fracture local to the operating site, elevation of segmentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count;
- mental illness;
- any patient unwilling to cooperate with postoperative instructions.

The contraindications of these devices are similar to those of other spinal interbody systems. This spinal instrumentation is not designed, or intended, or sold for uses other than those indicated.

Potential Adverse Events

In addition to the risks associated with spinal surgery without instrumentation, the following potential adverse events can occur (non-exhaustive list):

- delayed union;
- non-union (or pseudoarthrosis);
- bone fracture or stress shielding at, above, or below the level of surgery, microfracture, resorption, or damage of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site;
- retropulsed graft;
- loosening, or migration of implant, that may result in paralysis, neurological disorders or blood vessel erosion due to the proximity of the device;
- loss of neurological function, appearance of radiculopathy, dural tears, and/or development of pain, neurovascular compromise including paralysis, temporary or permanent retrograde ejaculation in males, or other types of serious injury, cerebral spinal fluid leakage;
- infection;
- hemorrhage of blood vessels and/or hematomas;
- discitis, arachnoiditis, and/or other types of inflammation;
- subsidence of implant into adjacent vertebral bodies resulting in loss of height and possible impingement of neural structures;
- loss of proper spinal curvature, correction, height, and/or reduction;
- pain, discomfort, or other abnormal sensations due to the presence of the device;
- herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery;
- loss of or increase in spinal mobility or function;
- pressure on the surrounding tissues or organs;
- foreign body reaction due to the presence of implants, such as a mass, auto-immune disease and/or impaired healing;
- cessation of any potential growth of the operated portion of the spine,
- deep venous thrombosis, thrombophlebitis;
- development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.;

- bone graft donor site complication;
- □ urinary retention or loss of bladder control or other types of urological system compromise;
- scar formation possibly causing neurological compromise or compression around nerves and/or pain;
- change in mental status;
- incapacity to resume the activities of normal everyday life;
- death;

NOTE: Some potential adverse events may require an additional surgical revision procedure.

Packaging

Packages for each of the components should be intact upon receipt. Care should be taken not to damage the implants when opening the packaging. Damaged packages or products must not be used and should be returned to SpineVision®.

Cleaning and Sterilization

The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF instruments are supplied not sterile. They must be sterilized before use.

The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF systems consist of clean but not sterile implants which must be sterilized before use.

All packaging and labeling materials must be removed prior to sterilization. Decontamination and cleaning must be completed prior to sterilization. For details regarding the following sections, please refer to instructions for cleaning, sterilization, inspection and maintenance (lubrication) of SpineVision® medical devices provided with the sets.

For cleaning and decontamination of the non-sterile implants and instruments at the Healthcare facility:

It is recommended to use a neutral pH, enzymatic and alkaline (pH<11) cleaning agents solution and deionized or distilled warm (room temperature) water for soaking, cleaning and rinsing.

Disassemble necessary instruments. Articulated instruments must be opened.

Inspect each device and repeat the washing in the presence of any residual debris, paying close attention to threads, cannulas, hinges and hard to reach areas.

Note: Certain cleaning solutions such as those containing bleach or formalin may damage some devices, and should not be used, except in the case of aluminium instruments used on high risk patients.

Lubrication specification:

SpineVision® recommends the lubrication of instruments with moving or articulating parts. The instruments must be lubricated after pre-disinfection/decontamination and before steam sterilization.

Sterilization: Prior to use, the SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF instruments and non-sterile implants should be steam sterilized according to the following parameters:

METHOD	CYCLE	TEMPERATURE	MINI EXPOSURE TIME	MINI DRYING TIME
Steam	Prevacuum	134°C (273°F) See note 1	18 minutes	30 minutes
Steam	Prevacuum	132°C (270°F) See note 2	4 minutes	30 minutes

Validations were performed with dynamic-air-removal steam sterilization process.

Note 1: This sterilization cycle is not considered by the Food and Drug Administration to be a standard sterilization cycle. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the Food and Drug Administration for the selected sterilization cycle specifications (time and temperature).

Note 2: Remove all the silicon handles from the trays for this cycle and sterilize them separately in a double wrap. The instruments and implants must be completely dried before use.

Verification: Devices must always be verified before use: any devices presenting signs of weakness or surface scratches must not be used. (Return to the manufacturer or distributor any implants showing deterioration)

IMPORTANT NOTE: As the PEEK glass transition temperature is 145°C (293°F), it is important not to use a higher sterilization temperature than 134°C (273°F). Otherwise the implant could suffer deformation.

To minimize the risk of deformation due to too high temperatures, do not place any load on top of the implants during cleaning and sterilization.

Warnings and precautions

PREOPERATIVE

- The use of SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF is to be done using dedicated adequate instrumentation by a trained orthopedic surgeon or neurosurgeon familiar with these recommendations as well as with the operating techniques relating to these implants.
- The surgeon must be aware of not only the medical and surgical aspects of the implant, but also of the material limitations of PEEK implants. The patient must be instructed in the limitations of the implant, especially with regard to the weight bearing and other body stresses on the device before a solid fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device. The patient must be instructed that device failure could lead to the need for additional surgery to remove failed components.
- Only patients that meet the criteria described in the indications should be selected and patients corresponding to the aforementioned contra-indications should not be selected. The following factors may be of importance when selecting patients for internal stabilization devices :
 - foreign body sensitivity
 - certain degenerative diseases
 - senility, mental illness, alcoholism, or drug abuse
 - obesity
 - patient's occupation or activity level
 - smoking

INTRAOPERATIVE

- Correct handling of the implants is extremely important. The implants must be transported in a way not to incur any damage. The implants must not be scratched and must not collide during handling.
- Care must be taken when picking up the SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF implants, filling it with autogenous bone graft or bone substitute and inserting it between the vertebral bodies. Alterations may produce defects in surface finish and may become the focal point for eventual breakage of the implant or for decreased efficacy in withstanding post-operative migration (e.g. when teeth of implant are damaged).

- Correct selection of the implant size is extremely important. The chance of adequate stabilization is increased by the selection of the proper size of the implant regarding the anatomy and conditions of the patient. Proper selection of the implant can minimize, but not eliminate risks.
- The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF components have been designed for use with this system. Do not use the SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF with any components not specifically recommended, or with any components from another manufacturer.
- The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF cannot be inserted via an anterior approach.
- The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF cannot be used in cervical and thoracic spine.
- The implant should be placed on the peripheral cortical bone of the endplates to minimize the risk of post-operative subsidence.
- Surgical implants must never be reused. An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage

POSTOPERATIVE

- Careful patient monitoring for the first two to four months following the operation is very important while the fusion mass matures and becomes able to share load with the implant. Adequately instruct the patient in the appropriate postoperative care. The patient must be instructed to reduce stress on the implants in order to avoid clinical problems that may accompany failed fixation. The patient's ability and willingness to follow instructions are one of the most important aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal stabilization devices by migrating, loosening or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may occur if excessive loading is placed on the implant. An active, debilitated, or demented patient who is unable to use orthotic weight bearing devices may be particularly at risk.
- No internal stabilization device can withstand activity levels equal to those withstood by normal healthy segment. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.
- Implants may fail when subjected to the increased loading associated with delayed union or non-union. Internal stabilization devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur or is delayed, the device may break due to fatigue. Patients should be informed of the risks of implant failure.

Magnetic Resonance (MR) Safety: The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF have not been evaluated for safety and compatibility in the MR environment. The SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF has not been tested for heating or migration in the MR environment.


Product complaints

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identification, durability, reliability, safety, effectiveness and/or performance, should notify SpineVision®. Further, if any of the implanted SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately.

If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

In Europe	In USA
SpineVision 10 Rue de la Renaissance 92160 Antony France Tél.: +33 (0) 1 5333 2525 Fax: +33 (0) 1 5333 2539	SpineVision Inc. 185 Alewife Brook Parkway - Suite 402 Cambridge, MA 02138 - USA Phone: +1 888 770 4118 Fax: +1 888 770 4119
Email: corp.quality@spinevision.com	

	See package insert for labeling limitation
	Do not use if package is damaged
	Do not re-use
	Keep dry
	Keep away from sunlight
	Federal law (US) restricts this device to sale by or on the order of a physician.
	Consult instructions for use

IMPORTANT INFORMATION CONCERNING SpaceVision® PLIF, SpaceVision® TLIF, and SpaceVision® OLIF Systems INSTRUMENTATION

Outside U.S.A. – ENGLISH VERSION

SpineVision® Spacevision® PLIF: Posterior Lumbar Interbody Fusion Cage

SpineVision® Spacevision® TLIF: Transforaminal Lumbar Interbody Fusion Cage

SpineVision® Spacevision® OLIF: Oblique Lumbar Interbody Fusion Cage

Purpose

The Spacevision® PLIF system (Posterior Lumbar Interbody Fusion), the Spacevision® TLIF system (Transforaminal Lumbar Interbody Fusion) and Spacevision® OLIF (Oblique Lumbar Interbody Fusion Cage) are intended as an internal spacer between the vertebral bodies of the lumbar (L1-S1 for PLIF, L1-L5 for TLIF) spine used in surgical procedures to give a primary stability of the anterior column (with or without restoring disc height) until fusion occurs, via a posterior approach.

The Spacevision® PLIF/TLIF/OLIF cages must be combined with an appropriate posterior lumbar fixation device, e.g. the SpineVision® UNI-Thread™ or PLUS® or LUMIS™.

The system is also intended to enhance the development of a solid fusion.

According to the clinical evaluation report, it can be concluded that the benefice-risk ratio is positive.

Description

The Spacevision® PLIF system (Posterior Lumbar Interbody Fusion), the Spacevision® TLIF system (Transforaminal Lumbar Interbody Fusion) and Spacevision OLIF system (Oblique Lumbar Interbody Fusion Cage) consist of cages of varying heights, widths, depths, angulations and ancillary products for placement of the cages.

For the Spacevision® PLIF system, two cages must be used per segment for fusion in order to stabilize the segment concerned.

For the Spacevision® TLIF system and Spacevision® OLIF, one cage must be used per segment for fusion in order to stabilize the segment concerned.

The Spacevision® PLIF/TLIF/OLIF must be filled with bone graft or a bone graft substitute to enhance the bone fusion.

The choice of graft will always be the surgeon's responsibility.

All cages are machined from implant grade PEEK-OPTIMA® LT1.

Indications and Usage

Properly used, the Spacevision® PLIF system (Posterior Lumbar Interbody Fusion), the Spacevision® TLIF system (Transforaminal Lumbar Interbody Fusion) and Spacevision OLIF system (Oblique Lumbar Interbody Fusion Cage) are intended to aid in the development of a spinal fusion.) These systems are intervertebral body fusion devices indicated for use with autogenous bone graft in skeletally mature patients with Degenerative Disc Disease (DDD) at one or two continuous levels from L2-S1. The surgeon must be aware of not only the medical and surgical aspects of the implant, but also of the material limitations of PEEK-OPTIMA® (Polyetheretherketone) implants. Patients should have received at least 6 months of non-operative treatment prior to treatment with Spacevision® PLIF, Spacevision® TLIF or Spacevision® OLIF systems. The patient must be instructed in the limitations of the implant (see Warnings and Precautions), especially with regard to the weight bearing and other body stresses on the device before a solid fusion occurs. The patient must be instructed that non-compliance with postoperative instructions could lead to failure of the device. The patient must be instructed that device failure could lead to the need for additional surgery to remove failed components..

The Spacevision® PLIF/TLIF/ OLIF are indicated for degenerative disc disease, lumbar pseudarthrosis, and spondylolisthesis.

Packaging

Packages for each of the components should be intact upon receipt. Care should be taken not to damage the implants when opening the packaging. Damaged packages or products must not be used, but should be returned to SpineVision®.

Cleaning and Sterilization

The Spacevision® PLIF/TLIF/ OLIF instruments and implants are supplied not sterile. They must be sterilized before use. All packaging and labeling materials must be removed prior to sterilization. Decontamination and cleaning must be completed prior to sterilization.

For the initial and any subsequent cleaning, the recommended general protocol for the pre-disinfection, decontamination and cleaning of the implants and instruments is:

Disassemble necessary implants and instruments. Articulated instruments must be opened.

In a pre-disinfecting bath, immerse and soak (minimum of 15 minutes) implants and instruments in a quaternary ammonia detergent. Transfer implants and instruments to an ultrasound bath, immerse and soak (15 minutes) implants and instruments under ultrasound action, in a quaternary ammonia detergent. Detergents must be used according to manufacturer recommendations. It is recommended to use deionized or distilled warm (room temperature) water for soaking, cleaning and rinsing.

After soaking, scrub all surfaces with a soft appropriate brush, paying close attention to threads, cannulas, hinges and hard to reach areas.

Rinse immediately and thoroughly after washing.

Inspect each device and repeat the washing in the presence of any residual debris.

Immediately dry products.

Note: Certain cleaning solutions such as those containing bleach or formalin may damage some devices, and should not be used, except in the case of aluminium instruments used on high risk patients (see Note below).

Lubrication specification:

SpineVision® recommends the lubrication of instruments that are mechanical or contain articulating surfaces. The lubricant used for the lubrication process of surgical instruments must be available in spray and is recommended to bear the CE marking according to the Directive 93/42/EEC modified by 2007/47/EEC. The instruments must be lubricated after each pre-disinfection/decontamination and before each steam sterilization.

Sterilization: Prior to use, the SpaceVision® PLIF, SpaceVision® OLIF, and SpaceVision® TLIF instruments and non-sterile implants should be steam sterilized according to the following parameters:

METHOD	CYCLE	TEMPERATURE	MINI EXPOSURE TIME
Steam	Vacuum	134°C (273°F)	18 minutes

The instruments and implants must be completely dried before use.

Check for evidence of deterioration on the implants after sterilization (return to the manufacturer or distributor any implants showing deterioration).

IMPORTANT NOTE: As the PEEK Optima® LT1 glass transition temperature is 145°C (293°F), it is important not to use a higher sterilization temperature than 134°C (273°F). Otherwise the implant could suffer deformation.

To minimize the risk of deformation due to too high temperatures, do not place any load on top of the implants during cleaning and sterilization.

Note for European audiences only - New requirements related to high risk patients:

In the case of high risk patients, i.e. being suspected of contact with Non-Conventional Transmissible Agents (e.g. prions), decontamination should be performed according to requirements of the World Health Organization. In this case, SpineVision® instruments can be decontaminated with sodium hydroxide, with the strict exception of instruments made out of aluminum, where bleach should be used.

Postoperative mobilization

Careful patient monitoring for the first two to four months following the operation is very important while the fusion mass matures and becomes able to share load with the implant. The patient must be instructed to reduce stress on the implants in order to avoid clinical problems that may accompany failed fixation.

Contraindications

- Trauma, tumor.
- severe osteoporosis, calcium metabolism disorder
- pregnancy
- infection

- other relative contraindications include, but are not limited to, obesity, certain degenerative diseases and foreign body sensitivity. The patient's activity level, mental condition, or occupation may be factors that affect the surgical outcome.
- other conditions, including alcoholism and drug abuse may also place excessive stresses on the device.
- do not use this device in the presence of any neural or vascular deficit or compromising pathology that may be further injured by device application.

Warnings

- the use of Spacevision® PLIF/TLIF/OLIF is to be done by a trained orthopedic surgeon or a trained neurosurgeon familiar with these recommendations as well as with the operating techniques, relating to these implants and using dedicated adequate instrumentation.
- correct selection of the implant size is extremely important. The chance of adequate stabilization is increased by the selection of the proper size of the implant regarding the anatomy and conditions of the patient. Proper selection of the implant can minimize, but not eliminate risks.
- no internal stabilization device can withstand activity levels equal to those withstood by normal healthy segment. No implant can be expected to withstand unsupported weight-bearing stress indefinitely.
- the Spacevision® PLIF/TLIF/ OLIF components have been designed for use with this system. Do not use the Spacevision® PLIF/TLIF/ OLIF with any components not specifically recommended, or with any components from another manufacturer.
- the Spacevision® PLIF/TLIF/ OLIF can not be inserted via an anterior approach.
- the Spacevision® PLIF/TLIF/ OLIF can not be used in cervical and thoracic spine.
- implants may fail when subjected to the increased loading associated with delayed union or non-union. Internal stabilization devices are intended as load-sharing devices until normal healing occurs. If normal healing does not occur or is delayed, the device may break due to fatigue. Patients should be informed of the risks of implant failure.
- patient selection: the following factors may be of importance when selecting patients for internal stabilization devices :
 - foreign body sensitivity
 - certain degenerative diseases
 - senility, mental illness, alcoholism, or drug abuse
 - obesity
 - patient's occupation or activity level
 - smoking

General Precautions

- surgical implants must never be reused. An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage.
- correct handling of the implants is extremely important. The implants must be transported in a way not to incur any damage. The implants must not be scratched and must not collide during handling. Care must be taken when picking up the Spacevision® PLIF/TLIF/ OLIF, filling it with bone graft or bone substitute and inserting it between the vertebral bodies. Alterations may produce defects in surface finish and may become the focal point for eventual breakage of the implant or for decreased efficacy in withstanding post-operative migration (e.g. when teeth of implant are damaged).
- adequately instruct the patient in the appropriate postoperative care. The patient's ability and willingness to follow instructions are one of the most important aspects of successful bone healing. The patient must be instructed in the limitations of the implant, and that physical activity and full weight bearing may cause premature failure of internal stabilization devices by migrating, loosening or fracture. The patient must be instructed that an implant does not have the strength of normal healthy bone, and that device failure may

occur if excessive loading is placed on the implant. An active, debilitated, or demented patient who is unable to use orthotic weight bearing devices may be particularly at risk.

- proper placement of the implant. The implant should be placed on the peripheral cortical bone of the endplates to minimize the risk of post-operative subsidence.

Possible Adverse Effects

- delayed union
- non-union
- fracture, loosening, or migration of implant
- infection
- subsidence of implant into adjacent vertebral bodies resulting in loss of height and possible impingement of neural structures
- pain, discomfort, or other abnormal sensations due to the presence of the device
- cage migration that may result in paralysis or blood vessel erosion due to the proximity of the device. Such erosion may lead to hemorrhage which could result in death.
- foreign body reaction due to the presence of implants, such as a mass, auto-immune disease and/or impaired healing
- arrest of growth of fused segments of the spine
- deep vein thrombosis, thrombophlebitis and/or pulmonary embolism
- incapacity to resume the activities of normal everyday life
- death

Product complaints

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identification, durability, reliability, safety, effectiveness and/or performance, should notify SpineVision®. Further, if any of the implanted Spacevision® PLIF/TLIF/ OLIF components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately.

If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence (email: corp.complaint@spinevision.com). When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

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Email: corp.quality@spinevision.com

IMPORTANT INFORMATION CONCERNING Spacevision® PLIF System and Spacevision® TLIF System and Spacevision® OLIF System INSTRUMENTATION

SpineVision® Spacevision® PLIF: Cage postérieure intersomatique lombaire de fusion
SpineVision® Spacevision® TLIF : Cage transforaminale intersomatique lombaire de fusion
SpineVision® Spacevision® OLIF : Cage oblique intersomatique lombaire de fusion

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