Lumis UliS

INSTRUCTION FOR USE

ENGLISH (OUT OF US).....	2
UNITED STATES OF AMERICA.....	6
FRANÇAIS	11
DEUTSCH	16
ITALIANO	26
ESPAÑOL	27
PORTUGUESE.....	32

Introduction:

The Universal Lumbar Intuitive System (U.L.I.S.TM System), and Lumbar Universal Minimally Invasive System (LUMISTM System) instrumentations are designed for correction and surgical stabilization of the spine during development of solid bone fusion. It is recommended to remove the device as soon as effective solid bone fusion has been achieved.

According to the clinical evaluation report, it can be concluded that the benefice-risk ratio is positive.

Description:

The Universal Lumbar Intuitive System (U.L.I.S.TM System), and Lumbar Universal Minimally Invasive System (LUMISTM System) instrumentations are composed of pedicle screws and fixation rods. Their components can be rigidly assembled in a variety of constructs, each corresponding to the needs and anatomy of a specific patient.

These constructs are assembled using specific instruments.

The components of the U.L.I.S.TM system are made from ASTM F136 titanium alloy (Ti-6Al-4V ELI) complying with ASTM F136 (ISO 5832-3) or ASTM F1537 Cobalt Chromium.

The components of the LUMISTM system are made from ASTM F136 titanium alloy (Ti-6Al-4V ELI) complying with ASTM F136 (ISO 5832-3)

Implants must never be reused.

Components of the U.L.I.S.TM and LUMISTM systems must not be used with components derived from another manufacturer.

Indications, Contraindications and Potential adverse effects:

Indications:

When used for anterior screw fixation or as a posterior, non-pedicle system of the non-cervical spine, the U.L.I.S.TM and LUMISTM systems are indicated for:

- degenerative disc disease (discogenic back pain with degeneration of the disc confirmed by history and radiographic studies)
- spondylolisthesis
- fracture
- spinal stenosis
- tumors
- failed previous fusion (pseudoarthrosis)

The U.L.I.S.TM and LUMISTM systems are pedicle screw systems indicated for skeletally mature patients who:

- have severe spondylolisthesis (Grades 3 and 4) at the L5-S1 vertebra;
- receive fusions using autogenous bone graft only;
- have the device fixed or attached to the lumbar and sacral spine (L3 to sacrum); and
- have the device removed after the development of a solid fusion.

In addition, the U.L.I.S.TM and LUMISTM systems are pedicle screw systems intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar, and sacral spine :

- degenerative spondylolisthesis with objective evidence of neurologic impairment
- fracture
- spinal tumor
- failed previous fusion (pseudoarthrosis)

This device can only be implanted by a surgeon with a good working knowledge of the device, its applications, the instruments and the required surgical technique.

Contraindications:

Contraindications include, but are not limited to:

- Allergy to the implanted material, mainly to metal (e.g. cobalt, chromium, nickel, etc.)
- Any other medical or surgical condition likely to compromise the success of instrumented surgery, such as the presence of a malignant tumour or serious congenital abnormalities, raised erythrocyte sedimentation rate not explained by other diseases, high white blood cell count or a tendency to low white blood cell count.
- All cases not described in the indications.
- Localized infection of the operative site.
- All patients with insufficient tissue cover of the operative site.
- Local signs of inflammation.
- Fever or leukocytosis.
- Pathological obesity.
- Pregnancy.
- Mental illness.
- Rapidly evolving joint diseases, bone absorption, osteopenia and/or osteoporosis. Osteoporosis is a relative contraindication, as this medical condition can limit the expected correction gain and stability of mechanical fixation.
- All cases not requiring bone graft or bone fusion.
- All cases requiring a combination of different metals.
- All patients not agreeing to comply with postoperative instructions.
- The use of pedicular screws above T10.

The contraindications of these devices are similar to those of other spinal rod instrumentations. This spinal instrumentation is not designed, or intended or sold for uses other than those indicated.

Potential adverse effects:

In addition to the risks associated with spinal surgery without instrumentation, the following potential adverse events can occur (non-exhaustive list):

- Early or delayed hardware removal.
- Implant migration.
- Dismantling, deformity, sliding and/or rupture of one or all components or instruments.
- Foreign body reaction due to the presence of implants, such as a mass, auto-immune disease, metallosis and/or impaired healing.
- Subcutaneous pressure by the components, possibly causing alteration of the skin in places in which the tissue cover is insufficient. Skin complications, including perforation of the skin by the implant or the graft.
- Loss of curvature or correction of the spine, loss of height.
- Infection.
- Fracture of the vertebral body due to transfer of loads above and below the instrumented zone.
- Lack of bone consolidation (or non-union).
- Loss of neurological functions, appearance of radiculopathies, dural injuries and/or pain. Neurovascular insufficiency, including paralysis or other serious lesions. Cerebrospinal fluid leak.
- Gastrointestinal, urological and/or reproductive tract disorders, including sterility, impotence.
- Haemorrhage and/or haematoma.
- Arrest of growth of fused segments of the spine.
- Disk inflammation, arachnoiditis and/or other types of inflammation.
- Deep vein thrombosis, thrombophlebitis and/or pulmonary embolism.
- Bone graft donor site complications.
- Incapacity to resume the activities of normal everyday life.
- Soft tissue damage.
- Death.

- **NOTE:** Some potential adverse events may require an additional surgical revision procedure.

Placement of the device:

The U.L.I.S.[™] and LUMIS[™] systems are implants designed to be connected by rods 6mm in diameter. They are designed to be used with instruments specifically intended for these devices.

Connections for the U.L.I.S.[™] and LUMIS[™] systems must be tightened to a value of 8,5 N.m whatever the rod material is.

Packaging:

The packaging of each component must be intact on reception. The operating room personnel must thoroughly verify, before use, that all devices are complete and that no component presents any signs of damage. Damaged packaging and/or products must not be used, but must be returned to SpineVision[®].

All packaging and labelling materials must be removed prior to sterilization. Only sterile implants and instruments must be used in surgery.

Cleaning and Sterilization:

The U.L.I.S.[™] and LUMIS[™] systems instruments and implants are supplied not sterile. They must be sterilized before use. Surface decontamination and cleaning of reusable instruments must be completed prior to sterilization.

For the initial and any subsequent cleaning, the recommended general protocol for the pre-disinfection, decontamination and cleaning of **the implants and instruments** is:

Disassemble necessary implants and instruments. Articulated instruments must be opened.

In a pre-disinfecting bath, immerse and soak (minimum of 15 minutes) implants and instruments in a quaternary ammonia detergent. Transfer implants and instruments to an ultrasound bath, immerse and soak (15 minutes) implants and instruments under ultrasound action, in a quaternary ammonia detergent. Detergents must be used according to manufacturer recommendations. It is recommended to use deionized or distilled warm (room temperature) water for soaking, cleaning and rinsing.

After soaking, scrub all surfaces with a soft appropriate brush, paying close attention to threads, cannulas, hinges and hard to reach areas.

Rinse immediately and thoroughly after washing.

Inspect each device and repeat the washing in the presence of any residual debris.

Immediately dry products.

Note: Certain cleaning solutions such as those containing bleach or formalin may damage some devices, and should not be used, except in the case of aluminium instruments used on high risk patients (see Note below).

Lubrication specification:

SpineVision[®] recommends the lubrication of instruments that are mechanical or contain articulating surfaces. The lubricant used for the lubrication process of surgical instruments must be available in spray and is recommended to bear the CE marking according to the 93/42/EEC Directive modified by 2007/47/EEC Directive.. The instruments must be lubricated after each pre-disinfection/decontamination and before each steam sterilization.

Prior to use, the U.L.I.S.[™] and LUMIS[™] systems components should be steam sterilised according to the following parameters:

METHOD	CYCLE	TEMPERATURE	MINIMUM EXPOSURE TIME
Steam	Vacuum	134°C (273°F)	18 minutes

The instruments and implants must be completely dried before use.

Verification:

Devices must always be verified before use: any devices presenting signs of weakness or surface scratches must not be used.

Note:

New requirements related to high risk patients:

In the case of high risk patients, i.e. being suspected of contact with Non-Conventional Transmissible Agents (e.g. prions), decontamination should be performed according to requirements of the World Health Organization. In this case, SpineVision® instruments can be decontaminated with sodium hydroxide, with the strict exception of instruments made out of aluminium, where bleach should be used.

WARNINGS:

This instrumentation is not designed to be the only means of long-term support of the spine. The use of this product cannot be successful without a mechanically solid bone graft. In the absence of a solid bone graft, the implanted devices can become deformed, loose, dismantled and/or may break.

The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to severe spondylolisthesis (grades 3 and 4) of the L5-S1 vertebra, degenerative spondylolisthesis with objective evidence of neurologic impairment, fracture, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other conditions are unknown.

Compliance with preoperative and intraoperative procedures and recommendations, a good knowledge of surgical techniques, correct selection and positioning of implants as well as the quality of the reduction obtained are important factors determining success of the operation. Appropriate patient selection and patient cooperation also have a major influence on the results. High non-fusion rates have been demonstrated in smokers, obese subjects, alcoholics, patients with poor quality bone or muscle and/or suffering from paralysis. These patients must be informed of this risk and its consequences.

In the case of a major bone defect of the anterior vertebral column, the surgeon must consider the use of additional support devices.

Patient selection:

- Particular attention must be paid to patient selection criteria, correct placement of the implant and postoperative care in order to minimize the stresses to which the implants are submitted.
- Only patients satisfying the criteria described in the indications must be selected.
- Patients corresponding to criteria described in the contraindications must not be selected.

Implant Use:

- Adequate selection of the type, shape and size of the implant for each patient is crucial to the success of the operation. After implantation, implants are subjected to repeated stresses and their resistance is limited by adaptation of their dimensions to bone anatomy. When these dimensions are not adapted, repeated stresses can induce excessive loading of the hardware, resulting in deformity, rupture or dismantling of the device before bone consolidation has been achieved, which may cause damage or require premature removal of the instrumentation.
- Implants must be handled and stored very carefully. They must not be scratched or damaged. Implants and instruments must be protected during storage, especially from corrosive environments.

- As this is a mechanical assembly, the surgeon must be familiar with all components before using this instrumentation and must personally manipulate the components to ensure that all implants and instruments are available before starting surgery.
- The type of construct to be performed for each case must be determined before starting the surgical operation. An adequate range of implant sizes must be available at the time of the operation, including smaller sizes and larger sizes than those initially planned.
- All components must be sterilized before use. Supplementary sterile components must be available if required unexpectedly.
- Mixing of dissimilar metals can accelerate the corrosion process. The components of this system are not to be used with implants of other material in building a construct. Components of this system should NOT be used with components from any other system or manufacturer.
- U.L.I.S.[™] and UNI-Thread[®] ASTM F136 Titanium Alloy (Ti-6Al-4V ELI) and ASTM F1537 Cobalt Chromium rods can be used with the U.L.I.S.[™] system.
- LUMIS[™] and UNI-Thread[®] Titanium rods can be used with the LUMIS[™] system

PRECAUTIONS:

Precaution: The implantation of pedicle screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

SURGICAL IMPLANTS MUST NEVER BE REUSED:

An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage.

Preoperative:

- Implants must be handled and stored very carefully. They must not be scratched or damaged. Implants and instruments must be protected during storage, especially from corrosive environments.
- The surgeon must be familiar with all components before using this instrumentation and must personally manipulate the components to ensure that all implants and instruments are available before starting surgery.
- The type of construct to be performed for each case must be determined before starting the surgical operation. An adequate range of implant sizes must be available at the time of the operation, including smaller sizes and larger sizes than those initially planned.
- All components must be sterilized before use. Supplementary sterile components must be available if required unexpectedly.

Intraoperative:

- The surgeon must strictly comply with the instructions for use of the spinal instrumentation.
 - The surgeon must always exert extreme caution in relation to the spinal cord and nerve roots. This warning is particularly important during insertion of screws. All nerve lesions can induce loss of neurological function.
 - Rupture, sliding or incorrect use of the instruments or components of the implant can injure the patient or operating room personnel.
- Rods must not be repeatedly or excessively bent beyond that which is absolutely necessary. Very thoroughly verify that the surfaces of the implant are neither scratched nor damaged. If the rods must be cut to a certain length, they must be cut to obtain a flat, blunt surface, perpendicular to the axis of the rod. **Cut the rod outside of the surgical field.**
- Do not to use screws of inappropriate dimensions (length, diameter) due to the risk of damage to the nerve roots or causing haemorrhage and/or avulsion.
 - Bone graft must be performed to ensure satisfactory fusion. Autologous bone grafts are essentially used with this instrumentation.
 - Before closure, all tightening screws must be tightened according to the instructions.

Postoperative:

- The postoperative advice given by the surgeon to the patient and the patient's compliance with this advice are extremely important.
- The patient must be informed about the limits of use of the device. If premature and/or excessive loading of the spine occurs prior to complete bone consolidation, the patient must be informed that complications such as deformity, dismantling and/or rupture of the device can occur.
- The patient or the device must not be exposed to mechanical vibrations which could induce dismantling of the device. The patient must be informed about this risk and must be advised to limit his or her physical activities, particularly lifting and torsion movements, as well as any sports activities. The patient must be advised to refrain from smoking and drinking alcohol during the process of consolidation of the bone graft.
- Patients must be informed that they will be unable to laterally flex their spine at the level of fusion and must be trained to compensate for this permanent physical limitation of their body's movements.
- Persistent absence of bone consolidation that cannot be immobilized places repeated and excessive stresses on the implant. Due to the mechanism of fatigue, these stresses can finally result in dismantling, deformity or rupture of the device. It is important to immobilize the zone of fusion and confirm this consolidation by radiological examination. If an abnormal absence of consolidation is observed or if the components become dismantled, deformed and/or broken, the device must be removed immediately, before it causes a serious lesion.
- Spinal implants are internal fixation devices designed to stabilize the operative zone during the normal consolidation process. Once consolidation has been achieved, these devices no longer have any useful function and can be removed by the surgeon. If the device is not removed after having completed its role, one of the following complications can occur :
 - corrosion, with a tissue reaction or localized pain,
 - migration of the implant resulting in a lesion,
 - risk of additional lesions due to postoperative trauma,
 - deformity, dismantling and/or rupture may make hardware removal difficult or impossible,
 - pain or abnormal feelings due to the presence of the device,
 - high risk of infection and
 - osteolysis due to transfer of mechanical loads.

All explanted devices must be treated in order to prevent reuse in another surgical procedure.

Product Complaints:

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identity, reliability, durability, safety, effectiveness and/or performance should notify SpineVision®. Further, If any of the implanted spinal instrumentation components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately. If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

In Europe
SPINEVISION
Regulatory Affairs Dpt
10 rue de la Renaissance
Bâtiment E
92160 Antony
FRANCE
Tél. : +33 (0) 1 5333 2525
Fax: +33 (0) 1 5333 2539
Email: corp.quality@spinevision.com

IMPORTANT INFORMATION CONCERNING U.L.I.S.™ System, and LUMIS™ System INSTRUMENTATION

Introduction:

The Universal Lumbar Intuitive System (U.L.I.S.TM System), and Lumbar Universal Minimally Invasive System (LUMISTM System) instrumentations are designed for correction and surgical stabilization of the spine during development of solid bone fusion. It is recommended to remove the device as soon as effective solid bone fusion has been achieved.

Description:

The Universal Lumbar Intuitive System (U.L.I.S.TM System), and Lumbar Universal Minimally Invasive System (LUMISTM System) instrumentations are composed of pedicle screws and Titanium fusion rods from the UNI-Thread[®] System.

Their components can be rigidly assembled in a variety of constructs, each corresponding to the needs and anatomy of a specific patient.

These constructs are assembled using specific instruments.

The components of the U.L.I.S.TM system are made of titanium alloy (Ti-6Al-4V ELI) complying with ASTM F136 (ISO 5832-3) or ASTM F1537 Cobalt Chromium.

The components of the LUMISTM system are made of titanium alloy (Ti-6Al-4V ELI) complying with ASTM F136 (ISO 5832-3).

Implants must never be reused.

Components of the U.L.I.S.TM and LUMISTM systems must not be used with components derived from another manufacturer.

Indications, Contraindications and Potential adverse effects:

Indications:

When used for anterior screw fixation or as a posterior, non-pedicle system of the non-cervical spine, the U.L.I.S.TM and LUMISTM systems are indicated for:

- degenerative disc disease (discogenic back pain with degeneration of the disc confirmed by history and radiographic studies)
- spondylolisthesis
- fracture
- curvatures (i.e. scoliosis, kyphosis, and/or lordosis)spinal stenosis
- tumors
- failed previous fusion (pseudoarthrosis)

The U.L.I.S.TM and LUMISTM systems are pedicle screw systems indicated for skeletally mature patients who:

- have severe spondylolisthesis (Grades 3 and 4) at the L5-S1 vertebra;
- receive fusions using autogenous bone graft only;
- have the device fixed or attached to the lumbar and sacral spine (L3 to sacrum); and
- have the device removed after the development of a solid fusion.

In addition, the U.L.I.S.TM and LUMISTM systems are pedicle screw systems intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar, and sacral spine (T10-S1):

- degenerative spondylolisthesis with objective evidence of neurologic impairment
- fracture
- curvatures (i.e. scoliosis, kyphosis, and/or lordosis)
- spinal tumor
- failed previous fusion (pseudoarthrosis)

This device can only be implanted by a surgeon with a good working knowledge of the device, its applications, the instruments and the required surgical technique.

Contraindications:

Contraindications include, but are not limited to:

- Allergy to the implanted material, mainly to metal (e.g. cobalt, chromium, nickel, etc.)
- Any other medical or surgical condition likely to compromise the success of instrumented surgery, such as the presence of a malignant tumour or serious congenital abnormalities, raised erythrocyte sedimentation rate not explained by other diseases, high white blood cell count or a tendency to low white blood cell count.
- All cases not described in the indications.
- Localized infection of the operative site.
- All patients with insufficient tissue cover of the operative site.
- Local signs of inflammation.
- Fever or leukocytosis.
- Pathological obesity.
- Pregnancy.
- Mental illness.
- Rapidly evolving joint diseases, bone absorption, osteopenia and/or osteoporosis. Osteoporosis is a relative contraindication, as this medical condition can limit the expected correction gain and stability of mechanical fixation.
- All cases not requiring bone graft or bone fusion.
- All cases requiring a combination of different metals.
- All patients not agreeing to comply with postoperative instructions.

The contraindications of these devices are similar to those of other spinal rod instrumentations. This spinal instrumentation is not designed, or intended or sold for uses other than those indicated.

Potential adverse effects:

In addition to the risks associated with spinal surgery without instrumentation, the following potential adverse events can occur (non-exhaustive list):

- Early or delayed hardware removal.
- Implant migration.
- Dismantling, deformity, sliding and/or rupture of one or all components or instruments.
- Foreign body reaction due to the presence of implants, such as a mass, auto-immune disease, metallosis and/or impaired healing.
- Subcutaneous pressure by the components, possibly causing alteration of the skin in places in which the tissue cover is insufficient. Skin complications, including perforation of the skin by the implant or the graft.
- Loss of curvature or correction of the spine, loss of height.
- Infection.
- Fracture of the vertebral body due to transfer of loads above and below the instrumented zone.
- Lack of bone consolidation (or non-union).
- Loss of neurological functions, appearance of radiculopathies, dural injuries and/or pain. Neurovascular insufficiency, including paralysis or other serious lesions. Cerebrospinal fluid leak.
- Gastrointestinal, urological and/or reproductive tract disorders, including sterility, impotence.
- Haemorrhage and/or haematoma.
- Disk inflammation, arachnoiditis and/or other types of inflammation.
- Deep vein thrombosis, thrombophlebitis and/or pulmonary embolism.
- Bone graft donor site complications.
- Incapacity to resume the activities of normal everyday life.
- Soft tissue damage.
- Death.

NOTE: Some potential adverse events may require an additional surgical revision procedure.

Placement of the device:

The U.L.I.S.[™] and LUMIS[™] systems are implants designed to be connected by rods 6mm in diameter. They are designed to be used with instruments specifically intended for these devices.

Connections for the U.L.I.S.[™] and LUMIS[™] systems must be tightened to a value of 8,5 N.m whatever the rod material is.

Packaging:

The packaging of each component must be intact on reception. The operating room personnel must thoroughly verify, before use, that all devices are complete and that no component presents any signs of damage. Damaged packaging and/or products must not be used, but must be returned to SpineVision[®].

All packaging and labelling materials must be removed prior to sterilization. Only sterile implants and instruments must be used in surgery.

Cleaning and Decontamination:

The U.L.I.S.[™] and LUMIS[™] systems instruments and implants are supplied not sterile. They must be sterilized before use. Surface decontamination and cleaning of reusable instruments must be completed prior to sterilization.

For the initial and any subsequent cleaning, the recommended general protocol for the pre-disinfection, decontamination and cleaning of **the implants and instruments** is:

Disassemble necessary implants and instruments. Articulated instruments must be opened.

In a pre-disinfecting bath, immerse and soak (minimum of 15 minutes) implants and instruments in a quaternary ammonia detergent. Transfer implants and instruments to an ultrasound bath, immerse and soak (15 minutes) implants and instruments under ultrasound action, in a quaternary ammonia detergent. SpineVision recommends the use of Hexanios G+R for this process with Ultrasonic vibrations, deionized or distilled warm (room temperature) water for soaking, cleaning and rinsing.

After soaking, scrub all surfaces with a soft appropriate brush, paying close attention to threads, cannulas, hinges and hard to reach areas.

Rinse immediately and thoroughly after washing.

Inspect each device and repeat the washing in the presence of any residual debris.

Immediately dry products.

Note: Certain cleaning solutions such as those containing bleach or formalin may damage some devices, and should not be used, except in the case of aluminium instruments used on high risk patients (see Note below).

Lubrication specification:

SpineVision[®] recommends the lubrication of instruments that are mechanical or contain articulating surfaces. The lubricant used for the lubrication process of surgical instruments must be available in spray and is recommended to bear the CE marking according to the 93/42/EEC Directive modified by 2007/47/EEC Directive. The instruments must be lubricated after each pre-disinfection/decontamination and before each steam sterilization.

Steam Sterilisation

Prior to use, the U.L.I.S.[™] and LUMIS[™] systems components should be steam sterilised with according to the following parameters:

Method	Cycle	Temperature	Mini exposure time	Mini drying time
Steam	Prevacuum	134°C (273°F) See note 1	18 minutes	30 minutes
Steam	Prevacuum	132°C (270°F) See note 2	4 minutes	30 minutes

Validations were performed with dynamic-air-removal steam sterilisation process.

Note 1: This sterilization cycle is not considered by the Food and Drug Administration to be a standard sterilization cycle. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the Food and Drug Administration for the selected sterilization cycle specifications (time and temperature).

Note 2: Remove all the handles from the trays for this cycle and sterilize them separately in a double wrap.

The instruments and implants must be completely dried before use.

Verification:

Devices must always be verified before use: any devices presenting signs of weakness or surface scratches must not be used.

Note:

New requirements related to high risk patients:

In the case of high risk patients, i.e. being suspected of contact with Non-Conventional Transmissible Agents (e.g. prions), decontamination should be performed according to requirements of the World Health Organization. In this case, SpineVision[®] instruments can be decontaminated with sodium hydroxide, with the strict exception of instruments made out of aluminium, where bleach should be used.

WARNINGS:

This instrumentation is not designed to be the only means of long-term support of the spine. The use of this product cannot be successful without a mechanically solid bone graft. In the absence of a solid bone graft, the implanted devices can become deformed, loose, dismantled and/or may break.

The LUMIS[™] and U.L.I.S.[™] systems have not been evaluated for safety and compatibility in the MR environment. The LUMIS[™] and U.L.I.S.[™] systems have not been tested for heating or migration in the MR environment.

The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to severe spondylolisthesis (grades 3 and 4) of the L5-S1 vertebra, degenerative spondylolisthesis with objective evidence of neurologic impairment, fracture, spinal tumor, and failed previous fusion (pseudoarthrosis). The safety and effectiveness of these devices for any other conditions are unknown.

Compliance with preoperative and intraoperative procedures and recommendations, a good knowledge of surgical techniques, correct selection and positioning of implants as well as the quality of the reduction obtained are important factors determining success of the operation. Appropriate patient selection and patient cooperation also have a major influence on the results. High non-fusion rates have been demonstrated in smokers, obese subjects, alcoholics, patients with poor quality bone or muscle and/or suffering from paralysis. These patients must be informed of this risk and its consequences.

In the case of a major bone defect of the anterior vertebral column, the surgeon must consider the use of additional support devices.

Patient selection:

- Particular attention must be paid to patient selection criteria, correct placement of the implant and postoperative care in order to minimize the stresses to which the implants are submitted.
- Only patients satisfying the criteria described in the indications must be selected.
- Patients corresponding to criteria described in the contraindications must not be selected.

Implant Use:

- Adequate selection of the type, shape and size of the implant for each patient is crucial to the success of the operation. After implantation, implants are subjected to repeated stresses and their resistance is limited by adaptation of their dimensions to bone anatomy. When these dimensions are not adapted, repeated stresses can induce excessive loading of the hardware, resulting in deformity, rupture or dismantling of the device before bone consolidation has been achieved, which may cause damage or require premature removal of the instrumentation.
- Implants must be handled and stored very carefully. They must not be scratched or damaged. Implants and instruments must be protected during storage, especially from corrosive environments.
- As this is a mechanical assembly, the surgeon must be familiar with all components before using this instrumentation and must personally manipulate the components to ensure that all implants and instruments are available before starting surgery.
- The type of construct to be performed for each case must be determined before starting the surgical operation. An adequate range of implant sizes must be available at the time of the operation, including smaller sizes and larger sizes than those initially planned.
- All components must be sterilized before use. Supplementary sterile components must be available if required unexpectedly.
- Mixing of dissimilar metals can accelerate the corrosion process. The components of this system are not to be used with implants of other material in building a construct. Components of this system should NOT be used with components from any other system or manufacturer.
- SpineVision ASTM F136 Titanium Alloy (Ti-6Al-4V ELI) and ASTM F1537 Cobalt Chromium rods can be used with the U.L.I.S.TM system.
- LUMISTM and UNI-Thread[®] ASTM F136 Titanium Alloy (Ti-6Al-4V ELI) can be used with the LUMISTM system

PRECAUTIONS:

Precaution: The implantation of pedicle screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

SURGICAL IMPLANTS MUST NEVER BE REUSED;

An explanted implant should never be reimplanted. Even though a device appears undamaged, it may have small defects and internal stress patterns that may lead to early breakage.

Preoperative:

- Implants must be handled and stored very carefully. They must not be scratched or damaged. Implants and instruments must be protected during storage, especially from corrosive environments.
- The surgeon must be familiar with all components before using this instrumentation and must personally manipulate the components to ensure that all implants and instruments are available before starting surgery.
- The type of construct to be performed for each case must be determined before starting the surgical operation. An adequate range of implant sizes must be available at the time of the operation, including smaller sizes and larger sizes than those initially planned.
- All components must be sterilized before use. Supplementary sterile components must be available if required unexpectedly.

Intraoperative:

- The surgeon must strictly comply with the instructions for use of the spinal instrumentation.
- The surgeon must always exert extreme caution in relation to the spinal cord and nerve roots. This warning is particularly important during insertion of screws. All nerve lesions can induce loss of neurological function.
- Rupture, sliding or incorrect use of the instruments or components of the implant can injure the patient or operating room personnel.
- Rods must not be repeatedly or excessively bent beyond that which is absolutely necessary. Very thoroughly verify that the surfaces of the implant are neither scratched nor damaged. If the rods must be cut to a certain length, they must be cut to obtain a flat, blunt surface, perpendicular to the axis of the rod. **Cut the rod outside of the surgical field.**
- Do not to use screws of inappropriate dimensions (length, diameter) due to the risk of damage to the nerve roots or causing haemorrhage and/or avulsion.
- Bone graft must be performed to ensure satisfactory fusion. Autologous bone grafts are essentially used with this instrumentation.
- Before closure, all tightening screws must be tightened according to the instructions.

Postoperative:

- The postoperative advice given by the surgeon to the patient and the patient's compliance with this advice are extremely important.
- The patient must be informed about the limits of use of the device. If premature and/or excessive loading of the spine occurs prior to complete bone consolidation, the patient must be informed that complications such as deformity, dismantling and/or rupture of the device can occur.
- The patient or the device must not be exposed to mechanical vibrations which could induce dismantling of the device. The patient must be informed about this risk and must be advised to limit his or her physical activities, particularly lifting and torsion movements, as well as any sports activities. The patient must be advised to refrain from smoking and drinking alcohol during the process of consolidation of the bone graft.
- Patients must be informed that they will be unable to laterally flex their spine at the level of fusion and must be trained to compensate for this permanent physical limitation of their body's movements.
- Persistent absence of bone consolidation that cannot be immobilized places repeated and excessive stresses on the implant. Due to the mechanism of fatigue, these stresses can finally result in dismantling, deformity or rupture of the device. It is important to immobilize the zone of fusion and confirm this consolidation by radiological examination. If an abnormal absence of consolidation is observed or if the components become dismantled, deformed and/or broken, the device must be removed immediately, before it causes a serious lesion.
- Spinal implants are internal fixation devices designed to stabilize the operative zone during the normal consolidation process. Once consolidation has been achieved, these devices no longer have any useful function and can be removed by the surgeon. If the device is not removed after having completed its role, one of the following complications can occur :
 - corrosion, with a tissue reaction or localized pain,
 - migration of the implant resulting in a lesion,
 - risk of additional lesions due to postoperative trauma,
 - deformity, dismantling and/or rupture may make hardware removal difficult or impossible,
 - pain or abnormal feelings due to the presence of the device,
 - high risk of infection and
 - osteolysis due to transfer of mechanical loads.
 -

All explanted devices must be treated in order to prevent reuse in another surgical procedure.

Product Complaints:

Any health care professional (e.g., customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identity, reliability, durability, safety, effectiveness and/or performance should notify SpineVision®. Further, If any of the implanted spinal instrumentation components does not meet any of its performance specifications or otherwise does not perform as intended, or is suspected of doing so, SpineVision® should be notified immediately. If any SpineVision® product may have caused or contributed to the death or serious injury of a patient, SpineVision® should be notified as soon as possible by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report is requested from SpineVision®.

FOR ALL INFORMATION OR COMPLAINTS, CONTACT:

SpineVision Inc.

185 Alewife Brook Parkway
Suite 402
Cambridge, MA 02138 – USA
Phone: + 1 888 770 4118
Fax: +1 888 770 4119

Email: corp.quality@spinevision.com

IMPORTANT INFORMATION CONCERNING U.L.I.S.™ System, and LUMIS™ System INSTRUMENTATION

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